11. FAULT FINDING

Fault Description	Action
Scooter drives slowly, and	Batteries are low. Recharge the batteries.
indicator light flashes slowly	
The wheels pull to one	Check the tire's air pressure.
direction	
Noises come from the motor	Make sure no foreign objects or dirt is stuck in
and/or the wheels	the wheels or wheel hubs. If you do not find any
	obstructions, contact the service center.
The scooter does not turn on	Check that the key is inserted properly.
	Make sure that the batteries are inserted in
	the right position.
	Check the battery charge.
Batteries do not charge	Check the connection between the charger
	and the scooter charge outlet.
	Make sure that the batteries are inserted in
	the right position.
	Make sure the charge indicator light is lit
	during charging
The scooter turns on, but	Make sure the wheel release handle is in
does not drive	drive position (See Fig. 4).
	Verify that there are no objects interfering
	with the scooter.
	Check the battery charge status.

If the scooter does not operate properly and the indicator light is flashing quickly, count the number of flashes in each sequence in order to define the problem in accordance with the table below.

Number of Flashes	Description
1	Either the battery needs charging or the batteries
	are not connected properly.
2	There is no connection to the motor. Check the
	connections between the motor and the controller.
3	There is a short circuit between the motor and the
	battery.
6	Check if the charger is connected.
7	Problem with the operation lever. Make sure it is
	set in the middle before operating the scooter.
8	Problem with the controller. Make sure all the
	connections into the control panel are tight and
	secure.
9	There is no connection between the controller and
	the brake. Make sure all connections to the brake
	and controller are tight and secure.
10	Controller overload. Likely a result of a failure in the
	battery connection. Check connections to the
	battery and controller.