

# About the DST

Stairs, Slope and Parallel-bars

The DST provide a comprehensive, tailor made gait rehabilitation practice for all three essential functions for patient's mobility- walking, climbing stairs, and negotiating slopes. By a click of a button, stairs' height and slope's degree can be easily adjusted to meets patients' current abilities and needs.



#### **DST**8000

- 4 height adjustable stairs, from 0 to 16.5 cm.
- In horizontal position 1.5 meter parallel bars.



#### DST8000 Triple

- 4 height adjustable stairs, from 0 to 16.5 cm.
- Adjustable slope, from 0 to 26 degrees.
- 3 meters parallel bars in horizontal position.



### DST Triple Sense 💷 📈 ---- 😪

- 4 height adjustable stairs, from 0 to 16.5 cm.
- Adjustable slope, from 0 to 26 degrees.
- 3 meters parallel bars in horizontal position.Computerized- with
- full gait analysis on all surfaces, to diagnose, treat, and study patients' ambulatory capabilities.



TZORA ACTIVE SYSTEMS LTD Original factory: Kibbutz Tzora - Israel Main office: Simtat Ha-levanon, Giv'at Ko'ah, 7318000 Israel Phone: +972-35502020 WhatsApp +972-50-4474467 info@tzora.com US WebSite: tzora.com UK WebSite: tzora.co.uk IL WebSite : tzora.co.il Medical Trainers: tele-apt.com



 $\square$ 

- **DST**8000 Pro
- 4 height adjustable stairs, from 0 to 16.5 cm.
- In horizontal position 1.5 meter parallel bars.
- Computerized to monitor, record, and display patient's performance in past and current treatment sessions.



## **DST**8000 Triple Pro

- 4 height adjustable stairs, from 0 to 16.5 cm.
- Adjustable slope, from 0 to 26 degrees.
- 3 meters parallel bars in horizontal position.Computerized to
- monitor, record, and display patient's performance in past and current treatment sessions.

